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Alfa Step iQ Surface Profiler Manual

2011 / 3 / 30, rev. 2011 / 5 / 11

ATTENTION :Stylus element will be damaged if you move sample w/o CONFIRMING STYLUS UP. Repair will cost 1 million yen!!

1. Start up

- Confirm air suspension.
- Turn ON power switch under desk.
- Turn ON PC and LCD.
- Startup Alfa Step iQ
- Select User name and Enter password.
(Each laboratory has own user name. If not, you can select "Engineer")
- Select Setting / Stylus Force Calib. tab to confirm the stylus force.**
- Write your name, lab name, and **stylus force** on the log book.

2. Set your sample

- CONFIRM STYLUS is UP.** (By your eyes and monitor on LCD)
- Move in your sample under the stylus. Use XY knobs on the left side.

3 Adjustment of measuring point

- Continuously push button (or page down key) to move the stylus downward.
When you double click (or End key), the stylus move to touch the sample.
- Press (or page up key) to lift up the stylus.
- CONFIRM STYLUS is UP.** Adjust measurement position by XY knobs. (To maximize monitor, select "Maximize Video" by right-clicking on the monitor.)

4 Measurement

- Press button to load your special recipe.
(If you do not have, load default recipe .Do not overwrite default recipe.)

Default Recipe

Scan Length :500 μm	Scan Direction :
Scan Speed :20 $\mu\text{m} / \text{s}$	Required Stylus Force :GREEN
Sampling Rate :50 Hz	Range / Resol. :20 $\mu\text{m} / 1.19 \text{ pm}$
Adjustment :Center	

- Press button to modify the recipe. (To save it, press)
- Press (or F8 key) to start the measurement.

5 Data processing

5-1 Level compensation (mechanical)

•Adjust leveling knob. (Expert only.)

5-2 Level compensation (software)

•Press **Level** button, and select left and right positions. Then press **OK**.

5-3 Step Height

•Move L and R indicators. (e.g. L1 and L2 for averaging window). Step height is appeared at "StpHt".

5-4 Print

•Press **Print** button.

5-5 Save data

•Press **Data** and select "Save Raw (Current) Data". Select destination folder to save pro format (AlfaStep iQ format).

5-5 CSV Export

•Press **Data Preview** and select "File / Export results into a new text file" to save CSV data. Then press **Scan** to return scan mode.

6. Exchange sample

•Double-click **↑** (or Home key) to move upward the stylus.

CONFIRM STYLUS is UP. (By your eyes and monitor on LCD)

Take out your sample. Use XY knobs on the left side.

•To measure the next sample, go to "2. Set your sample" on this manual

7. Shut down

•Press **Exit** to shut down.

•Turn OFF LCD and power switch under the desk.

•Close air valves.

•Write end time on the log book

8. Miscellaneous

Prohibited for general users:

Stylus force adjustment and stylus replacement

If necessary, please ask _____ (ext. 3136).